


Overview

Powerful new *in vitro* assays provide a translational method to study biologics or small molecule modulators of immune responses. We've developed an optimized panel of T cell assays available to help clients better understand the complex biology of immuno-oncology.



DISCOVERY

Click to learn more


Complex Biology *In Vitro* Assays: Immuno-Oncology Cytokine Response Assay

TCR-Dependent Stimulation of Th-1, Th-2 and Th-17 Cytokine Biosynthesis in Isolated Human PBMCs

The immune system has a very important role in tumor progression; hence, it is crucial to characterize the possible impact of a drug candidates (antibody or small molecule therapeutic) on the immune system. Our [immune cell activation assays](#) assess the impact of candidate compounds on multiple super antigen stimulated T-cell receptor (TCR) engagement pathways. The [immuno-oncology assays](#) are optimized for multiplex cytokine analysis and identified cytokines can be future studied *in vivo* across an [immunology platform](#).

Assay Principle

Freshly isolated human PBMCs from healthy donors are seeded in the absence or presence of various stimuli (plate-bound anti-CD3, soluble anti-CD28, or PHA) for the indicated time (48 hours) in the presence or absence of test compounds. After 2 days, the cell culture supernatant is removed for multiplex cytokine analysis.



Need a custom version of this assay?

Visit criver.com/ds-vitro-assay

EVERY STEP OF THE WAY


Need a custom version
of this assay?
Visit criver.com/ds-vitro-assay

Assay Setup

T cell-activated protocol has been developed for optimum analysis of Th-1/Th-2/Th-17 cytokine response.

T Cell Activation Assay	
Donor	PBMCs
Seeding density	2×10 ⁵ c/w
Trigger (positive control)	anti-CD3 alone, anti-CD3 + anti-CD28 or PHA
Incubation	48 hours post-trigger
Readout	Cytokine release

Representative dose response data shown below from one donor with multiple stimulants measured at 48 hours post-stimulation. The cytokines have been chosen to represent Th-1 (IFN-γ and IL-2), Th-2 (IL-5, IL-10 and IL-13) and Th-17 (IL-17) responses. Where both anti-CD3 and anti-CD28 stimulation were performed, a dose response to anti-CD3 was tested in the presence of a fixed concentration of anti-CD28 (10 ug/mL), except the graph for IL-2, in which case 5 ug/mL of anti-CD28 was used.



Need a custom version
of this assay?
Visit criver.com/ds-vitro-assay

Summary

Type 1 T helper (Th-1) cells produce proinflammatory cytokines interferon-gamma, interleukin (IL)-2, and tumor necrosis factor (TNF)-beta, which activate macrophages and are responsible for cell-mediated immunity and phagocyte-dependent protective responses. By contrast, type 2 T helper (Th-2) cells produce IL-4, IL-5, IL-10, and IL-13, which are responsible for strong antibody production, eosinophil activation, and inhibition of several macrophage functions, thus providing phagocyte-independent protective responses. Th-2 responses will counteract the Th-1 mediated action.

Th-17 cells are another T helper subset that are characterized as preferential producers of interleukin-17A (IL-17A), IL-17F, IL-21, and IL-22. Th-17 cells and their effector cytokines mediate host defensive mechanisms to various infections, especially extracellular bacterial infections, and are involved in the pathogenesis of many [autoimmune diseases](#). The optimal scenario would therefore be that there should be a well-balanced Th-1, Th-2 and Th-17 response, suited to the specific immune challenge.

Here we show robust dose responses from multiple donors evaluating key cytokines that are involved in these T helper cell cytokine responses, Th-1 (IFN- γ and IL-2), Th-2 (IL-5, IL-10 and IL-13) and Th-17 (IL-17) using three different stimuli exploring [T cell activation pathways](#).

Using these [immuno-oncology assays](#) with specified stimuli of interest, changes in immunoreactivity of molecules can be identified in isolated primary human PBMCs from our in-house donor panel. A single cytokine or a panel can be chosen from the examples shown to evaluate immunoreactivity of molecules.

Complementary Immuno-Oncology Assays

[T Cell Proliferation CTG Assay](#)

[T Cell Exhaustion Assay](#)

[T Cell-Mediated Chemotaxis Assay](#)

[3D Spheroid T Cell Cytotoxicity Assay](#)