

Profiling and Heightened Characterization of O-linked Glycans




Elaine Sun¹, Mario DiPaola¹, Marshall Bern², Andrew Hanneman¹
¹Charles River Laboratories, Woburn, MA; ²Protein Metrics Incorporated, San Carlos, CA


1 Key Points

- Biopharma requires robust methods for routine profiling and detailed structural characterization of O-linked glycans
- Monitoring the degradation reactions associated with maintaining an intact reducing end may dictate the most suitable approaches
- Advantageous if: the same method used to profile oligosaccharides includes monosaccharide O-glycans in the same chromatogram
- The exact methods for pursuing detailed structural analysis (accurate mass, MS/MS, and MSⁿ) is an important consideration, especially for complex glycoproteins such as blood derived protein products

2 O-glycan Diversity




Profiling and Characterization of O-glycans from Diverse/Complex O-glycoproteins




- ### Non-CHO Glycoproteins
- Emerging Expression Systems:** O-glycosylation in *Spodoptera frugiperda* (Sf9) and *Trichoplusia ni* (Hi-5) insect cell lines is complex and include abundant hexuronic acid (Sf9 and Hi-5) and O-linked phosphocholine (Sf9)
Glycobiology vol. 22 no. 11 pp. 1593 Abstract #217
Stefan Gaunitz, Chunsheng Jin, Anki Nilsson, Jining Liu, Niclas G. Karlsson and Jan Holgersson
 - Glycoengineering:** e.g. GlycoSwitch®
Jacobs PP, Geysens S, Vervecken W, Contreras R, Callewaert N. Engineering complex-type N-glycosylation in *Pichia pastoris* using GlycoSwitch technology. Nat Protoc. 2009; 4:58–70.
Hamilton SR, Cook WJ, Gomathinayagam S, Burnina I, Bukowski J, Hopkins D, Schwartz S, Du M, Sharkey NJ, Bobrowicz P, Wildt S, Li H, Stadheim TA, Nett JH. Production of sialylated O-linked glycans in *Pichia pastoris*. Glycobiology. 2013;23:1192–1203

3 MS-based Methods used to Profile and Characterize O-linked Glycans


- Analyze O-glycans while they still reside on the protein or peptide backbone (may need to de-N-glycosylate, de-sialylate, etc.)**
- Release by non-reductive β-elimination for subsequent chromophore or fluorophore labeling (important to monitor degradation)**



3. Release using “classical” reductive β-elimination (minimal degradation “peeling” product)




B) Native Alditols (PGC-LC-CAD-MS/MS)




4 LC-MS with Charged Aerosol Detection (CAD)

- Configuration**
- Release optimization**



3. Applications

- A) Biopharma samples by PGC-CAD-MS (CHO derived)**
- B) Blood derived O-glycosylated products (MSⁿ characterization of permethylated O-glycans)**



5 Conclusions

1. Currently, we prefer traditional reductive beta elimination due to the diversity of structures potentially presented in a contract research organization (CRO) environment (potential degradation reactions cannot be tolerated)
2. We closely monitor peeling reactions during method development
3. Place monosaccharides within the same chromatogram as oligosaccharides
4. CAD detection may be promising for routine batch-to-batch comparability
5. Prefer MS/MS for structural confirmation (vs. accurate mass)
6. Permethylated and MSⁿ are highly useful for detailed or *de novo* structure analysis

Acknowledgements: We would like to thank the entire CRL Woburn Mass Spectrometry Team: Kurt Morgenstern, Le Meng, Art Cansizoglu, and Michael Chi; as well as St. John Skilton, Rose Lawler and Ilker Sen from PMI.