

Meeting the Challenges of Antibody Drug Conjugate Characterization by LC-MS/(MS)

Jing Li¹, Andrew Hanneman, Mario DiPaola
Charles River, 8 Henshaw Street, Woburn, MA 01801
Mario.DiPaola@CRL.com 617-234-0001

¹Currently at Alnylam Pharmaceuticals Inc.

LakePharma's Protein Engineering Symposium

10/2017

www.criver.com

For access to additional scientific posters, please visit www.criver.com/posters.

Meeting the Challenges of Antibody Drug Conjugate Characterization by LC-MS/(MS)

Jing Li¹, Andrew Hanneman, Mario DiPaola



Charles River, 8 Henshaw Street, Woburn, MA 01801 Mario.DiPaola@CRL.com 617-234-0001

ABSTRACT


Antibody-drug conjugates (ADCs) are potent biopharmaceuticals comprised of a cytotoxic agent attached to an antibody or antibody fragment via a chemical linker. As complex three part (bio)molecules, ADCs require analytical considerations intersecting the sometimes disparate modes of small and large molecule investigations, and mass spectrometry plays a key role in their characterization. We describe LC-MS/(MS) approaches addressing various ADC attributes, including drug-antibody ratio (DAR) determination, mapping of drug conjugation sites, and accounting for various drug/linker chemistries including associated impurities. Selected MS-based methodologies are designed based on the class of ADC molecule, the conjugation chemistry, and any requirements for fine-scale structural understanding.

Analysis of ADCs with attachments through cysteine, lysine, or other sites may require significantly different strategies. The application of native vs. denaturing LC-MS, including instrument optimization in response to various linker chemistries, is discussed. For characterization of site-to-site drug conjugation, sequence-dependent selection of appropriate endopeptidases and application of appropriate MS/MS fragmentation approaches is considered. Impurity characterization is carried out by accurate mass and MS/MS analysis to a significant degree, and the limitations of MS-based methodologies are kept in mind with regard to the known challenges of small molecule structural characterization.

Results and Discussion

Characterization of Reduced Antibody Drug Conjugates with LC-MS

Figure 1. (A) UV 280 chromatograms of reduced unconjugated antibody and drug conjugated antibody separated on a hydrophobic interaction column (HIC) using an MS-friendly mobile phase. Deconvoluted mass spectra of (B) light and (C) heavy chains obtained with a high resolution mass QTOF spectrometer.



ADCs were deglycosylated, reduced, and the light and heavy chains, with varied degrees of lysine-conjugation, were well separated on a HIC column with the eluate analyzed by QTOF MS.

Drug-Antibody Ratio Determination via Intact LC-MS

Figure 2. DAR Analysis at the intact mass level under denaturing conditions: Acetonitrile /Water with TFA.





Figure 3. Analysis of Intact ADCs under denaturing vs. non-denaturing LC-MS conditions: (A) denaturing, (B) non-denaturing for cysteine conjugations.





CONCLUSIONS

Mass Spectrometry is a key tool for ADC characterization, including DAR determination, drug linked site mapping and investigations of drug-linker chemistry including related impurities.

DAR and Conjugation Site Analysis via LC-MS Peptide Mapping

Figure 4. Endopeptidase digestion followed by QTOF LC-MS/MS. (A) Identification of drug-linked peptides, (B) DAR tabulation.



ACKNOWLEDGMENTS

We would like to acknowledge all the members of Charles River for their support.
1: Currently at Alnylam Pharmaceuticals Inc.
The work described herein is intended for research purposes only; not intended for any human or animal therapeutic or diagnostic use, unless otherwise stated. The trademarks or copyrights mentioned herein are the property of their respective owners.

Impurity Characterization

Figure 5. Characterization of impurities in a peptide/protein conjugate. (A) Full-view C4 LC-MS chromatogram, (B) Zoom-in of impurities, (C) Identification by accurate mass, (D) by MS/MS, and (E) isotope patterns.

