

Adenovirus-mediated CRISPR/Cas9 gene editing in patient primary fibroblasts to validate potential drug targets in lung fibrosis

Anne-Marie Zuurmond^{1#}, Liesbeth Veugen¹, Lieke Geerts¹, Barbara Killian², Isabella Gazzoli¹, Dennis Gravekamp¹, Monique Hazenoot¹, Miranda van der Ham¹, Britt Sotthewes¹, Brigitta Witte¹, Krista Ouwehand¹, David Fischer³, Jeroen DeGroot¹

¹Charles River Leiden, Darwinweg 24, 2333 CR Leiden, The Netherlands; ²Charles River Wilmington, 251 Ballardvale St, Wilmington, MA 01887, USA; ³Charles River Cambridge, Chesterford Research Park, Saffron Walden, Essex CB10 1XL, United Kingdom
#Corresponding author: Anne-Marie.Zuurmond@crl.com



1 INTRODUCTION

Targeted gene editing has become a key tool in deciphering the molecular mechanisms underlying disease processes and for target validation in the drug discovery process. For research into new drug targets and compounds for treatment of pulmonary fibrosis, we previously developed a cell-based assay with primary lung fibroblasts derived from idiopathic pulmonary fibrosis (IPF) patients. In this in vitro model, the transition of fibroblast into myofibroblast (FMT) is triggered by TGF β . Myofibroblasts are thought to play a major role in fibrosis through excessive deposition of extracellular matrix. A key feature of myofibroblasts is the production of α -smooth muscle actin (α SMA).




Figure 1: Cell-based phenotypic assay for pulmonary fibrosis

In this study, we have evaluated the feasibility to validate potential drug targets in fibrosis by CRISPR/Cas9-mediated targeted gene knockout in this assay. For delivery of Cas9 and gRNA the adenoviral system was chosen that is known for its ability to transduce a wide range of cells with high efficiency, including human primary lung fibroblasts that are difficult to transfect. The ACTA2 gene encoding α SMA was selected as target to demonstrate feasibility of CRISPR/Cas9-mediated knockout in lung fibroblasts. The setup was validated using a small gRNA library with putative drug targets for FMT/IPF emerging from literature.

2 EXPERIMENTAL DESIGN

An adenoviral system for CRISPR/Cas9 delivery was constructed with Cas9 and the gRNA being expressed from two different adenoviruses (AdV). Four gRNAs targeting ACTA2 and ≥ 2 gRNAs per putative drug target were designed and cloned into the adapter vector of the adenoviral system. High titer stocks of adenoviruses expressing either Cas9 or ACTA2 gRNAs were produced and used in the transduction of primary fibroblasts derived from four IPF patients. CRISPR/Cas9-mediated knockout of ACTA2 was measured by high content analysis (HCA).




Figure 2: Timeline for a cell-based phenotypic assay for lung fibrosis. Fibroblasts derived from IPF patients were transduced with AdV-Cas9 and AdV-gRNA one day after seeding. The fibrotic process was triggered by adding 2 ng/ μ l TGF β 4 days after adenoviral transduction. The fibroblast to myofibroblast transition (FMT) was quantified by measuring α SMA production by high content analysis (HCA) 2 days after adding the TGF β trigger.

3 RESULTS




Figure 3: Impaired α SMA production in FMT assay due to CRISPR/Cas9-mediated ACTA2 knockout. IPF patient-derived fibroblasts from four donors (IPF03, -05, -07 and -08) were transduced at a MOI of 24 each for AdV-Cas9 and AdV-gRNA_ACTA2 (ACTA2_g). TGF β was added 4 days after transduction and α SMA production was quantified by HCA 2 days after TGF β trigger (day 7). A non-targeting control adenovirus (AcGFP_g1) was used to normalize the data as adenoviral transduction is known to stimulate α SMA production. All four gRNAs were effective in ACTA2 gene knockout in two independent experiments, however gene editing efficiency was donor-dependent. Typical high content images are shown.

| Gene | Full name | Rationale |
|--------|--|---|
| ACTA2 | Actin alpha 2, smooth muscle | Involved in vascular contractility, is used as a FMT marker |
| IRF5 | Interferon regulatory factor 5 | Transcription factor modulating cell growth and differentiation |
| ITGA5 | Integrin subunit alpha 5 | Involved in the formation of the fibronectin receptor, cell surface adhesion and signaling |
| ITGB7 | Integrin subunit beta 7 | Adhesion receptor involved in cellular signaling with the extracellular matrix and MAPK-ERK pathway |
| ITGB1 | Integrin subunit beta 1 | Involved in the formation of the fibronectin receptor and tissue repair processes |
| COL3A1 | Collagen type III alpha 1 chain | Fibrillary collagen found in extensible connective tissue, is used as a FMT marker |
| TGFBR3 | Transforming growth factor beta receptor 3 | Membrane proteoglycan often functioning as a co-receptor with other TGF β receptor superfamily members |
| SMAD4 | SMAD family member 4 | Signal transduction protein involved in TGF β signaling regulating transcription processes |
| MAPK8 | Mitogen-activated protein kinase 8 | Regulation of kinases signaling processes involved with cell proliferation and differentiation associated with renal fibrosis and fatty liver disease |
| CXCL10 | C-X-C motif chemokine ligand 10 | Ligand for the receptor CXCR3 and inflammation marker, involved with adhesion molecule expression |
| S100A8 | S100 calcium binding protein A8 | Involved in the regulation cell cycle progression and differentiation associated with cystic fibrosis |

Table 1: Selection of putative drug targets for modulation of α SMA protein levels in the FMT process. Candidate genes are selected that are presumed to be involved in IPF based on available literature. The ACTA2 gene is selected as an assay control. For each target gene ≥ 2 gRNAs were designed and cloned into the adapter vector of the adenoviral system.




Figure 5: Representative high content images of CRISPR/Cas9-mediated ACTA2, IRF5, ITGA5, COL3A1 and SMAD4 knockout using the top-performing gRNAs in donor IPF07. Relative percentage α SMA / nuclei levels from high content analysis is given normalized to the non-targeting control adenovirus (AcGFP_g1).




Figure 4: CRISPR/Cas9-mediated downregulation of α SMA protein in FMT assay was observed using gRNA directed against IRF5, ITGA5, COL3A1 and SMAD4. The relative effect size for modulation of each of these targets was donor- and gRNA-dependent and varied between the two experiments. Percentage α SMA / nuclei levels is given normalized to the non-targeting control adenovirus (AcGFP_g1). Non-triggered, TGF β -triggered, crude lysate/TGF β -treated and ACTA2 gRNA-exposed fibroblasts were taken along as controls. Dark green: $\leq 50\%$; light green: 51-75%; grey: $\geq 76\%$.

4 CONCLUSION

We conclude that adenoviral transduction of primary fibroblasts for targeting gene editing using CRISPR/Cas9 technology is sufficiently efficient to assess the effects of a gene knockout in a cell-based assay such as the FMT assay without the need for selection or cloning. This was validated not only by targeting the ACTA2 gene encoding for α SMA protein directly, but also using a small library of gRNAs targeting genes that are associated with the fibrosis process and/or IPF phenotypes. Using this approach IRF5, ITGA5, COL3A1 and SMAD4 are identified as modulators of myofibroblast transition in this assay. Such effective gene targeting in a functional assay is a valuable asset for target validation in drug discovery.